Basic Shade Composition

Shade	A1	A2	A3	A3.5	A4	B1	B2	B3	B4
Opaque	OPA0	A20	A30	A3.50	A40	B10	B20	B30	B40
Cervical	-	-	AC1	-	AC2	-	-	BC1	BC2
Opaque Dentin	ODA1	ODA2	ODA3	ODA3.5	ODA4	ODB1	ODB2	ODB3	ODB4
Body	A1B	A2B	A3B	A3.5B	A4B	B1B	B2B	B3B	B4B
Incisal	58		59		60	57	58	59	60

Shade	C1	C2	C3	C4	D2	D3	D4
Opaque	C10	C20	C30	C40	D20	D30	D40
Cervical	-	-	CC1	-	CC2	DC1	DC2
Opaque Dentin	ODC1	ODC2	ODC3	ODC4	ODD2	ODD3	ODD4
Body	C1B	C2B	C3B	C4B	D2B	D3B	D4B
Incisal	58		59		60	60	59

EN Light-curing Crowns and Bridge, Inlays and Veneer Composite

READ CAREFULLY BEFORE USE

CERAMAGE is a light-curing Crown and Bridge composite featuring esthetic shades similar to natural dentition, optimal viscosity for detailed reproduction and easy-to-use property.

COMPOSITION

PRE-OPAQUE: UDMA, Aluminium Silicate, 2-HEMA, Pigments and others
OPAQUE: UDMA, Aluminium Silicate, 2-HEMA, Glass Powder, Pigments and others
C&B COMPOSITE RESIN: UDMA, UDA, Zirconium Silicate, Pigments and others
FLOWABLE COMPOSITE RESIN: UDMA, UDA, Zirconium Silicate, Pigments and others
MODELLING LIQUID: UDMA, Trimethylolpropane trimethacrylate, others

INDICATIONS

- > Fabrication of permanent and temporary indirect restorations
 - a) Crowns and Bridges for anterior and posterior restorations
 - b) Jacket crowns for anterior and posterior restorations
 - c) Inlays and Onlays
 - d) Laminate veneers
 - e) Implant superstructures
- > Intraoral and extraoral repairs of indirect restorations

PRECAUTIONS

1. Do not use product on patients with known allergies to this material and/or methacrylate monomer.
2. Operators with known allergies to this material and/or methacrylate monomer must not use this product.
3. If any inflammation or other allergic reactions occur on either the patient or operator, immediately discontinue use and seek medical advice.
4. Do not use this product on patients with malocclusion and bruxism (clenching, grinding and tapping).
5. CERAMAGE SPACER includes natural rubber. Operators with history of developing anaphylactic reaction to natural rubber must not use CERAMAGE SPACER.
6. Do not handle unpolymerized material with bare hands. Use protective plastic gloves and glasses to avoid sensitization to this product. Avoid contact of this material with skin and eyes. In case of accidental contact with skin, immediately blot with alcohol moistened cotton ball and rinse with plenty of water. In case of contact with eyes, immediately flush the eyes with plenty of water and seek medical advice.
7. Use local dust extractor, dust protective mask, etc while grinding this product to avoid harmful influence of the dust on the human body.
8. Give patients proper instructions for daily cleaning to avoid getting stained or plaque adhesion.
9. This product is intended for use by dental professionals only.

DIRECTIONS FOR USE

* Light-curing time described in the DIRECTIONS FOR USE below is based on the curing time with SHOFU's Solidilite V, Halogen light-curing unit for lab use.

A. Fabrication of anterior and posterior crowns and bridges/implant superstructures

A-1. Production of a metal frame
The construction of metal framework is modelled anatomically and made according to the generally accepted guidelines. Apply retention beads in a size of 150 µm before casting. To obtain a natural tooth shape and sufficient stability, ensure a minimum thickness of 0.6 mm of the veneer. Adjust the metal framework after casting.

> Occlusal contact of posterior restoration must be supported by metal frame.

A-2. Surface treatment of metal frame
Sandblast the metal frame. Then, clean with a steam cleaner or ultrasonic cleaner and dry. Apply M.L. Primer to the surface of the metal framework, where CERAMAGE is applied, and dry.

A-3. Application and light-curing of PRE-OPAQUE/SHOFU Universal Opaque (PRE-OPAQUE)
Apply a coat of PRE-OPAQUE with a brush and light-cure for 1 minute. Eject the material completely flows into the undercuts of the retention beads.

> For a bridge case with pontic, apply FLOWABLE COMPOSITE RESIN/ CERAMAGE UP to make up for the cutback area. Then, light-cure for 3 minutes.

A-4. Application and light-curing of OPAQUE/SHOFU Universal Opaque (OPAQUE)

Apply the material with a small brush and light-cure for 3 minutes. If necessary, repeat this procedure to mask the metal color completely.

> Use separate brushes for PRE-OPAQUE and OPAQUE and clean them with alcohol (or CERAMAGE CLEANER) immediately after use.

A-5. Build up and light-curing of C&B COMPOSITE RESIN (Opaque Dentin, Cervical and Body)
Apply Opaque Dentin or Cervical from cervical area to the center of the crown and light-cure for 1 minute (preliminary curing). When build up Body and light-cure for 1 minute (preliminary curing). When additional build-up is required, apply MODELLING LIQUID thinly before build-up to facilitate the procedure.

A-6. Build-up and light-curing of C&B COMPOSITE RESIN (Incisal, Translucent and Concentrate)
Create crown outline with Incisal. When characterizing partially, use Translucent or Concentrate. For detailed characterization, use FLOWABLE COMPOSITE RESIN/ CERAMAGE UP to minimize entrapping air bubbles and facilitate operation. After final build-up, light-cure for 3 minutes. Applying CERAMAGE OXY-BARRIER before final cure is recommended in order to minimize uncured material on the surface.

A-7. Contouring
After light-curing, contour the restoration following the conventional method. Use Dura-Green Stones or fine diamond burs for trimming and contouring. It is not advisable to use cutters or course diamond burs.

A-8. Finishing and polishing
After finishing with CompoMaster Coarse, use DURA-POLISH for pre polishing and DURA-POLISH DIA for high luster polishing.

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B-1. Preparation on a working model
After creating the working model in conventional method, apply CERAMAGE SPACER to the die. Do not apply to the margin area. Then, apply CERAMAGE SEP to the margins etc and dry.

B-2. Application and light-curing of OPAQUE/SHOFU Universal Opaque (OPAQUE)
Apply the material with a small brush and light-cure for 3 minutes. If necessary, repeat this procedure to mask the abutment color completely.

B-3. Build up C&B COMPOSITE RESIN light-cure, contour, finish and polish, following the procedures of A-4 to A-8 steps described above.

C. Fabrication of inlays and onlays/laminate veneers

C-1. Production on a working model
Create a working model following the conventional method. In the case, where cavity has undercuts, apply CERAMAGE SPACER. Then, apply CERAMAGE SEP to the inner surfaces around the cavity or on the area to be veneered and dry.

C-2. Build-up and light-curing of C&B COMPOSITE RESIN (Body and Incisal)
Build up Body from the cavity floor and light-cure for 1 minute (preliminary curing). When additional build-up is required, apply MODELLING LIQUID thinly before build-up to facilitate the procedure. Then, overlay Incisal. When characterizing partially, use Translucent or Concentrate. After final build-up, light-cure (preliminary curing) for 1 minute. After final build-up, light-cure for 3 minutes. Applying CERAMAGE OXY-BARRIER before final cure is recommended in order to minimize uncured material on the surface.

C-3. Contour, finish and polish following the procedures of A-7 and A-8 steps described above.

Light-curing times

Procedure	Halogen lamp (Solidilite V)
Curing of PRE-OPAQUE	1 min.
Curing of OPAQUE	3 min.
Preliminary curing	1 min.
Curing of pontic	3 min.
Final light-curing	3 min.

> The table shows light-curing time with Solidilite V. When using other curing unit, follow the curing unit manufacturer's instructions for use to light-cure CERAMAGE properly.

NOTE ON USE

1. Average light-curing time for intraoral repair with CERAMAGE

Dental light-curing unit	Halogen	LED
Light-curing time	40 sec.	20 sec.

2. When using a hand-held light-curing unit, operators must use protective glasses or light-shielding plate to avoid looking directly at the curing light. Operators must also protect patients' eyes from direct exposure to the curing light.

3. Follow the instructions for use of each dental material, instrument or equipment to be used in conjunction with this product.

4. Tightly close the cap immediately after each use.

5. The paste dispensed onto the dish or paper pad should be shielded from ambient light by light-shielding cover to avoid unintended polymerization of the material.

6. The pastes should not be mixed each other or with other material in order to minimize entrapping air bubbles and degrading material properties.

7. FLOWABLE COMPOSITE RESIN is intended to compensate air bubbles or fine adjustment of the contour. Avoid building up full contour or applying to large area.

8. Do not use this product in places with intense light such as close to a window or under the work light to avoid unintended polymerization of the paste.

9. Do not use this product for any purposes other than specifically outlined in the INDICATIONS in these instructions for use.

10. Use this product within the expiration date indicated on the package and container. (Example ☞ YYYY-MM-DD → Year-Month-Date of the expiration date)

STORAGE

Store in dark place at room temperature (1-30 °C/34-86 °F). Avoid high humidity and keep away from direct sunlight and any source of ignition.

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

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