



TRENDGOLD
Dental Legierungen

Lot 1 PF 800

CE0123

Binder Dental GmbH, Werner-v-Siemens-Str. 41, 49124 Georgsmarienhütte
Telefon +49 5401 48080 Telefax +49 5401 480828 www.trendgold.com



Arbeitstemperatur in °C: 780
Schmelzintervall in °C: 690 - 800

Massentanteile in %:

Au	Ag	Pd	Pt	Zn	Cu	Ir	In	Rh
70,5	6,5	-	1,4	3,0	13,0	x	5,5	-